


AIR-FLOW MASTER®

→ THE WORLD'S FIRST SUBGINGIVAL
PROPHYLAXIS UNIT



reddot design award
winner 2008

INNOVATION
AWARD
2007




EMS 
ELECTRO MEDICAL SYSTEMS

**TOUCH
'N'
FLOW**



> Original AIR-FLOW® method developed for application in deep periodontal pockets



→ 2 APPLICATION SYSTEMS IN ONE

- > PERIO-FLOW® handpiece and PERIO-FLOW® nozzle with powder chamber for AIR-FLOW® powder PERIO
- > AIR-FLOW® handpiece and powder chamber for AIR-FLOW® powder CLASSIC or AIR-FLOW® powder SOFT

PERFECT FOR SUBGINGIVAL AND SUPRAGINGIVAL PROPHYLAXIS

ORIGINAL AIR-FLOW® METHOD → SUBGINGIVAL

PERIO-FLOW® HANDPIECE →

- > Innovative design, specially developed for subgingival application, down to deep periodontal pockets
- > Ergonomic handpiece for the new AIR-FLOW® powder PERIO
- > Sterilizable at 134 °C

PERIO-FLOW® NOZZLE →

- > 3 nozzle outlets to optimize gentle turbulence of air-powder mixture in the subgingival region, horizontal
- > 1 nozzle outlet for water, vertical
- > Patented single-use nozzle
- > Quick connection onto the new PERIO-FLOW® handpiece

AIR-FLOW® POWDER PERIO →

- > Grain size ~ 25 µm (d50)
- > Ultrafine grain matched to the new PERIO-FLOW® handpiece, with PERIO-FLOW® nozzle
- > Does not abrade dentin

SUBGINGIVAL APPLICATION OF THE ORIGINAL AIR-FLOW® METHOD REDUCES PERIODONTAL POCKET DEPTH, REMOVES BIOFILM, PREVENTS PERIIMPLANTITIS.



EMS-SWISSQUALITY.COM

For more information > welcome@ems-ch.com