

# Power Tower Silence 200 – network-compatible and convenient



COMPRESSED AIR  
SUCTION  
IMAGING  
DENTAL CARE  
HYGIENE

## Compressed air and suction in a single housing – now with continuous monitoring

- Easy to use with great visualization thanks to a 4.3" colour touch screen
- Continuous monitoring of compressed air and suction systems via PC
- Connection to surgery network
- Volume-optimized thanks to intelligent regulation of compressor generators
- Signals from the amalgam separator are displayed

# Power Tower Silence – Performance made-to-measure

The Power Tower Silence combines compressed air and suction in a single housing. The compact, quiet supply unit can be installed flexibly in the surgery on a floor space of only 0.4 m<sup>2</sup>. This saves valuable space compared with the installation of single appliances – a big advantage with present-day rental costs. Metal construction and housing meet high aesthetic demands and can be integrated into all interior design styles.

At the top the compressed air,  
at the bottom the suction

- Cabinet module PTS 200 with noise-absorbing housing
- VS 1200 S Combination Suction Unit  
factory-fitted into Cabinet Module
- CA 4 Amalgam Separator with 16 l/min throughflow  
and automatic centrifuge braking every 15 minutes
- Membrane-Drying Unit
- 2x Duo Compressor Aggregates



The 4,3" colour touch screen enables the monitoring directly at the Power Tower. Thanks to its network compatibility it can be monitored via the PC directly at the reception. All important parameters are constantly under control and should any failure take place they are reported immediately.

