

Dental Unit WaterLine—Colony result register

Please record your DUWL results here & monitor for re-test intervals with the MHPC TEST samplers

1. For initial test results below 200CFM a three monthly testing period is recommended, provided the effective DUWL maintenance procedures are adhered to. Maintenance should also be in conjunction with Adec ICX water treatment tablets on each water bottle fill.
2. If any results record over “200CFM”, a DUWL shock treatment will be required. After a shock treatment & 3 consecutive months of passing test results (below 200CFM each month) then you can revert back to the once every three month testing period.
3. These procedures are strongly recommended for use with Adec DUWL—with separate clean water bottle systems. These procedures may be applicable to other brands that have the similar manufacturer recommendations & DUWL system.

Records

Chair S/N I.D	Date of Test	CFU Count	Shock Treatment?	Retest date on...

Ivoclar Vivadent Ltd

12 Omega Street, Rosedale, Auckland 0632
Tel 09 914 9999 | Fax 09 914 9990 | Freephone 0508 IVOCLAR (486 252)
Email info.nz@ivoclarvivadent.com | www.ivoclarvivadent.co.nz