

Simply a good feeling – System Hygiene from Dürr Dental

Suction systems



The original since 1965

Orotol plus | Disinfectant for suction systems



Foam-free liquid concentrate for the simultaneous disinfection, cleaning, deodorisation and care of dental suction systems and amalgam separators

Spectrum of activity

Bactericidal, yeasticidal, limited virucidal (enveloped viruses incl. HBV, HCV, HIV and coronaviruses), as well as adenoviruses and noroviruses

**Active agents
Units**

Quaternary ammonium compounds
1 litre (50 doses) (771812)
2.5 litres (125 doses) (716448)

Ideal addition

MD 555 cleaner

Effectively dissolves prophy powder

MD 555 cleaner | Special detergent for suction systems



Foam-free intensive cleaning concentrate for all types of suction systems and amalgam separators

Ingredients

Organic and inorganic acids, foam-free surfactants

Units

1 litre (20 doses) (771813)
2.5 litres (50 doses) (716441)

Ideal addition

Orotol plus

OroCup | Care system

User-friendly, closed care system for disinfection and cleaning of suction systems and amalgam separators



Volume

2 litre care system (704711)

Ideal addition

Orotol plus, MD 555 cleaner

The video shows you how to correctly use Orotol plus and MD 555 cleaner:
www.duerrdental.com/how-to-orotol



i Why is a special detergent like MD 555 cleaner needed?

While Orotol plus removes organic deposits from the system, the special detergent MD 555 cleaner eliminates inorganic residues. This includes limescale, for example, but also prophy powder, which is being used more frequently these days. Without using MD 555 cleaner, the powder can quickly accumulate and block the system.



Reduced alcohol content

FD 366 sensitive | Cleaning & disinfection of sensitive surfaces

Alcohol-reduced, pH-neutral solution for the quick cleaning & disinfection of sensitive surfaces of invasive and non-invasive medical devices and fixtures

Spectrum of activity

Bactericidal, tuberculocidal, yeasticidal, limited virucidal (enveloped viruses incl. HBV, HCV, HIV and coronaviruses), noroviruses

Active agents

17% Alcohol, quaternary ammonium compounds

Units

2.5 litres (716447)
Dispenser of 100 wipes (20 x 30 cm) (737693)
Refills of 4x 100 wipes (737694)
Package with 50 wipes (17.5 x 25 cm) (759506)

Ideal addition

Dry wipes, FD 360



How to use FD 366

FD 322 | Quick-acting disinfection

Fast-acting, ready-to-use solution for the disinfection of alcohol-resistant surfaces of invasive and non-invasive medical devices and fixtures

Spectrum of activity

Bactericidal, tuberculocidal, fungicidal, limited virucidal (enveloped viruses incl. HBV, HCV, HIV and coronaviruses), limited virucidal PLUS (incl. rotaviruses, adenoviruses, noroviruses), polyomaviruses SV40

Active agents

Alcohols

Units

2.5 litres (727734)

Ideal addition

Dry wipes



FD 360 Cleaning and care of vinyl upholstery

Cleaning fluid for deep cleaning and care of vinyl upholstery

Ingredients

Special surfactants, silicon compounds and care formula based on avocado oil Set with 250 ml spray bottle, 1 FD polishing wipe and 5 special sponges (716434)

Units

Ideal addition

FD 366 sensitive



How to use FD 360