PEOPLE HAVE PRIORITY

Sterilization, Hygiene & Maintenance

Simple, fast, efficient.



Oil service in record time.
No waiting. No stress.

assistina TWIN

Oil service without waiting

Green light for added efficiency: With the new Assistina TWIN waiting time during the oil service is a thing of the past. It is the first reprocessing device to be equipped with two process chambers that can be used alternately, thus allowing non-stop loading. Relaxed maintenance more effective than ever!

> Wow, it's quick! The dual chamber system means I only need a fraction of the time for the maintenance of the instruments. It is extremely quiet and looks great too.
>
> Sandra Schnäp,
> dental assistant



Running smoothly

1



Connect the instrument in chamber 1.

2



Close chamber 1 and start process.

3



Connect the second instrument and close chamber 2 as soon as the process in chamber 1 is complete.

4




Remove instrument from chamber 1 and repeat steps 1 and 2.



Oil service in record time

The Assistina TWIN performs your instrument maintenance quicker than ever before. Through the innovative oil nebulization technology, W&H has been able to reduce the maintenance process to a record time of just 10 seconds. For lightning quick and incomparably thorough results.



Groundbreaking technology
Thorough, comprehensive oil service of gearing components (orange) and cleaning of spray channels (green) thanks to innovative oil nebulization technology.



An innovative system
The oil is first nebulized before being blasted through the instrument at high pressure. This removes contamination from the inside of the instrument even more effectively. The fine mist this produces coats even the most remote parts of the gears without any need for the gearing components to be set in motion first.
All aerosols are eliminated from the process chambers via active suction and an integrated HEPA filter.

Practical for every practice

The Assistina TWIN adapts to the specific requirements of your practice. It can process up to 360 instruments per hour with the innovative processing procedure with the short processing time and also the ergonomic design with a wide range of adaptors.





Adaptors for all situations
The two chambers can be equipped with corresponding adaptors to suit your requirements.

Now I have the certainty that my instruments can be reprocessed between treatments without additional waiting times for me or my patients. It's a win-win situation for my practice.

Dr. Christian Müller,
dentist with his own practice

Maximum oil service with minimal consumption

Thanks to the smart process monitoring, the Assistina TWIN always performs the maintenance thoroughly and with minimal consumption. The Assistina TWIN Care Set does not need to be replaced until after approximately 2,800 instruments. The maintenance of your instruments can be performed not only reliably and quickly, but also more cost-effectively than ever before.



TWIN Care Set
Comprises a cartridge of cleaning solution, a cartridge of Service Oil F1 and a HEPA filter. All the consumables can be replaced with ease and without the need for any tools.



Assistina TWIN Care Set

contains as much as 20 oil spray cans*





*400 ml per oil spray cans



Hygiene of the highest standard



Instruments
Straight, contra-angle and turbine handpieces, air-driven scalers, air motors

Thermo washer disinfectors
Cleaning
Disinfection

Assistina TWIN
Cleaning of spray channels
Oil service

Seal²
Packaging

Multidem
Water treatment

Lisa / Lina
Sterilization

LisaSafe
Labelling

Lisa mobile app
Traceability



Technical data

Scan the QR code with your mobile phone for more information on Assistina TWIN

CE

Type:	**MB-302 Assistina TWIN**
Process:	Cleaning of spray channels, oil service
Suitable for:	Straight, contra-angle and turbine handpieces, air-driven scalers, air motors
Supply voltage:	100 – 240 V AC
Maximum power consumption:	18 VA
Air consumption:	Approx. 40 Nl/min.
Operating pressure:	5 – 10 bar (controlled via an integrated, automatic pressure controller)
Dimensions (W x D x H):	297 x 435 x 190 mm
Weight:	3.5 kg
Capacity:	200 ml W&H Active Fluid cleaning solution 200 ml W&H Service Oil F1

PEOPLE HAVE PRIORITY

Oil service in record time

Thorough, efficient, safe



Perfect oil service in just 10 seconds

The new Assistina TWIN reprocesses your instruments quicker and more thoroughly than ever before. Thanks to the innovative oil nebulization technology, the unit can process up to 360 instruments per hour. The ideal, practical solution for every practice.

assistina TWIN

Pictures are for illustrative purposes only. Additional equipment and accessories shown are not included as standard.

W&H Dentalwerk Bürmoos GmbH
Ignaz-Glaser-Straße 53, Postfach 1
5111 Bürmoos, **Austria**
t + 43 6274 6236-0
f + 43 6274 6236-55
office@wh.com
wh.com

20458 AEN Rev. 000 / 15.02.2017
Subject to alterations